


# Automated Frozen Sample Aliquotter

## No-Thaw Automated Extraction of Multiple Frozen Aliquots from One Frozen Sample

**CRYOTRACT**

Dale Larson<sup>1,2</sup>, John Slusarz<sup>2</sup>, Steve Bellio<sup>1</sup>, Linda Maloney<sup>1</sup>, Helena Judge Ellis<sup>3</sup>, Nader Rifai<sup>4</sup>, Gary Bradwin<sup>4</sup>, Jorge de Dios<sup>5</sup>

<sup>1</sup>Draper Laboratory; <sup>2</sup>Harvard Medical School; <sup>3</sup>Brown University; <sup>4</sup>Children's Hospital, Boston; <sup>5</sup>CryoXtract Instruments, LLC



### introduction

Biotechnology tools hold potential for exponential progress in molecular medicine and biomarker discovery. Yet, quality results depend on quality materials, and optimal preservation and sampling are crucial to ensure banked sample integrity and promote ideal conditions for analysis.

Current methods and tools force biobanking into sample storage and processing trade-offs which can impact operations and costs significantly; sometimes affecting sample quality. For example, repeated freeze/thaw cycling may degrade critical biological molecules, and minimizing sampling frequency is recommended to reduce potential sample molecular damage.

### the technology

The Automated Frozen Sample Aliquotter enables the hands-free extraction of multiple frozen aliquots from one single vial of frozen serum or plasma sample without thawing it. A proprietary robotic probe cores a frozen sample under ultra-cold conditions and deposits the still-frozen cores into separate cryotubes for downstream analyses. The extracted cores remain frozen through the process. Any sample remaining may be returned to the freezer still frozen.




Robotic Coring Probe

### case studies

#### productivity gains


##### scenario

- Clinical biobank
- Samples frozen in 1.8mL/2.0mL tubes
- Request for 2 aliquots X150 samples = 300 aliquots



##### tasks


- Process 1 source/2 destinations
- Sub-aliquot remaining sample material
- Distribute samples
- Return unused sample to freezer



### frozen sample aliquotter

Conceived and designed to serve the needs of modern biobanking, the Automated Frozen Sample Aliquotter helps optimize frozen sample integrity and sampling efficiency:


- Extract multiple frozen cores from one single cryotube of frozen plasma or serum
- Deliver hands-free, automated operation (e.g. source vial selection, de-capping/capping, coring, dispensing, cleaning) after sample loading
- Achieve volumetric consistency on every extracted core (never <100µL from 1.8mL cryotube)
- Maintain ultra-cold conditions pre-, during, post-coring
- Avoid carryover between samples



### demonstrated capabilities

Independent evaluation of a proof-of-concept instrument at the R.I. BioBank successfully demonstrated that it can extract multiple frozen, uniformly-sized, consistently homogeneous frozen cores from one frozen plasma sample stored in a 1.8mL cryotube. The extracted cores give reproducible results with very low variability when analyzed for common analytes.

#### Consistently Homogenous Cores with Low Variability



NOTE: All results were normalized using the assay results from the controls

### collection & storage optimization

##### scenario



##### approach

- Freeze sample in multiple small aliquots (pre-aliquot)

#### Current Process

Pre-aliquot



#### Introduce aFSA

Freeze in larger volume



- 18 racks "savings" (51% "savings")
- Up-front labor
- Sample tracking
- Cost of Plastics
- QA

Integrity • extend sample usable life • standardize sample aliquotting • enhance lab productivity