

# Five-Year, Multi-Institution Validation Studies of Archived Amplified Genomic DNA, Methylated DNA, and cDNA: Enhancing Tissue Archives With Nucleic Acid Archives

J. Langmore, V. Makarov, E. Kamberov, T. Kurihara, E. Bruening, T. Tesmer, J. M'Mwirichia, Rubicon Genomics, 4370 Varsity Drive, Ann Arbor, MI 48108, (734) 368-1705 [langmore@rubicongenomics.com](mailto:langmore@rubicongenomics.com) [www.rubicongenomics.com](http://www.rubicongenomics.com)


## ABSTRACT

Epidemiologic, pharmacologic, and diagnostic studies link clinical with molecular data via analysis of the nucleic acids, proteins, and metabolites in patient samples. Although clinical data are immortal, clinical specimens are not. While the number and power of molecular studies grow exponentially, the numbers of clinical specimens grow very slowly due to the expense of recruiting and following up patients as well as depletion and degradation of existing samples. Typically 1-4 mL of serum or plasma, less than 200 mL of urine and between 1-10 mm of diseased tissue are archived from each patient, and have a shelf-life of less than 20 years. The quantity and quality of investigator research is severely compromised by the great difficulty in accessing the limited quantities of patient material, high cost and irreproducibility of cutting, staining, characterizing, and macro dissecting tissue samples, or of purifying molecules from individual sections or fluid aliquots.


Greater utilization of archived samples for genetic, epigenetic, and expression studies can be achieved by amplifying, aliquoting, and storing nucleic acids from small portion of each patient sample. Isolation, 10,000-fold amplification, and QC of the DNA, methylated DNA or cDNA costs only \$1 per microgram of amplified product—a small fraction of the cost of recruiting the patient and clinical data. Aliquots of amplified DNA are stored without degradation and can be shared with thousands of researchers for over 100 years, thus extending and increasing the long-term value of the patient sample and data. After DNA amplification, individual biofluid aliquots or microtome sections can be used repeatedly and reproducibly by different researchers doing complementary genetic, epigenetic, and expression studies using qPCR, microarrays and next-generation sequencing. These amplified nucleic acid archives are needed to meet the increasing demands of these high throughput analytical platforms.

Published and unpublished data from over five years of experience in pharmaceutical, diagnostic, university, and government laboratories using more than 30,000 frozen and fixed tissue, serum, plasma, and urine samples will be presented to document the precision, accuracy, clinical significance, and economics of analyzing amplified DNA, meDNA, and cDNA from cancer and normal patients.

## Archived Biospecimens Can Be Amplified as DNA, Methylated DNA, and cDNA to Create Inexhaustible Archives of Nucleic Acids



## Amplification of Frozen or Fixed Tissue, or Biofluid is a Robust, Automatable 3-Step Process



1. Amplification works equally well for intact or highly degraded nucleic acids.
2. Nucleic acid can be reamplified a billion-fold without loss of representation.
3. Extremely low background allows single cell and single chromosome amplification.

## Practical Advantages of Nucleic Acid Archives

- Only a small fraction of the biospecimen needs to be extracted and amplified to produce DNA and RNA for hundreds of investigators.
- Microtomy, dissection, DNA/RNA extraction, and quality control are needed only once per biosample.
- Archived aliquots of amplified DNA, meDNA, and cDNA can be retrieved and supplied to investigators within one day and for as little as \$1/microgram.
- Exactly the same DNA or cDNA sample can be studied by different investigators at different times.
- Archived amplified DNA, meDNA, and cDNA stable for decades, without change in signal or background.

## GenomePlex Whole Genome Amplification

Research kits sold by Sigma-Aldrich  
Services and diagnostic kits sold by Rubicon

### Archiving of GenomePlex WGA DNA

- Library synthesis takes ~ 2 h; 1,000X amplification takes 1 h.
- Amplified DNA performs as well as gDNA for genotyping, CGH, and sequencing.
- Low background allows genetic studies of single cells and chromosomes.
- GenomePlex amplifies highly degraded DNA (e.g., FFPE, serum/plasma, urine) for all applications except genotyping with long STRs.
- Excellent array profiling possible from 20 ng unfixed DNA or 200 ng fixed DNA.
- WGA DNA is stable if stored at -20 C, and can be reamplified >1E09 times without significant loss of representation—even for very GC-rich sequences.
- A 20 ng aliquot of biospecimen DNA can be amplified to >20 mg to supply DNA to thousands of investigators
- GenomePlex in current use to discover disease genes and biomarkers, to manufacture MDx products, and for diagnosis of human disease.

## MethylPlex Whole Methylome Amplification

Diagnostic Partnerships and Beta Kits Available from Rubicon

### Archiving of MethylPlex WMA DNA

- Library synthesis takes ~ 3 h; 1,000X amplification takes 1 h.
- NO bisulfite conversion; no methylation-specific PCR assays to design; no custom arrays necessary.
- Same clinical sensitivity and specificity as best bisulfite-based tests.
- 2 mL of plasma
- Amplified MethylPlex DNA used to assay for DNA methylation using any copy number assay, such as qPCR, promoter and CpG island arrays, or Next-Generation sequencing platforms.
- Low background allows methylation studies of as little as 10 cells, and as little as 0.01% methylation
- MethylPlex works on highly degraded DNA (e.g., FFPE, serum/plasma, urine) for all methylation assays.
- Excellent array methylation profiling from 20 ng unfixed DNA or 200 ng fixed DNA.
- MethylPlex WMA DNA is stable if stored at -20 C, and can be reamplified >1E06 times without significant loss of sensitivity.
- A 20 ng aliquot of biospecimen DNA can be amplified to >20 mg of MethylPlex DNA to enable thousands of investigators to study methylation in any part of the genome
- MethylPlex has been used to discover biomarkers, and develop patient tests.


## TransPlex Whole Transcriptome Amplification

Research kits sold by Sigma-Aldrich  
Services and diagnostic kits sold by Rubicon

### Archiving of TransPlex WTA cDNA


- Library synthesis takes ~ 2 h; 1,000X amplification takes 1 h.
- Much faster and simpler than IVT.
- Amplified cDNA performs as well as unamplified cDNA for qPCR assays or profiling of gene or exon expression.
- No 3' bias or exon bias, as found after IVT.
- TransPlex amplifies highly degraded RNA (e.g., FFPE, serum/plasma, urine) for all applications.
- Excellent array profiling possible from 20 ng unfixed RNA or 200 ng fixed RNA.
- WTA cDNA is stable if stored at -20 C, and can be reamplified >1E09 times without significant loss of representation.
- A 20 ng aliquot of biospecimen RNA can be amplified to >20 mg to supply DNA to thousands of investigators
- TransPlex has been used to discover disease genes and biomarkers.

## CONCORDANCE BETWEEN Q-PCR EXPRESSION OF GENES BEFORE AND AFTER TRANSPLEX AMPLIFICATION (Chinniyayan lab)




## NSABP GENE EXPRESSION PROFILING OF B-27.2 PRE-TREATMENT FFPE CORE BIOPSY SECTIONS TO PREDICT RISK OF BREAST CANCER RECURRENCE (Dr. Soon Paik, recent update to NCI)

Study of 326 core-biopsy cases show that expression of ESR1 has much better prognostic value than pathological complete response, node status, or TRT. STUDY DESIGN— •100 ng total RNA as starting material •RNA amplified to 20 µg using Rubicon TransPlex custom amplification service •Hybridization to Affymetrix GeneChip U133 2.0 plus •PAM and SUPERPC used for prediction of ER, pCR, and outcome. Microarray data was 96% concordant with IHC. ESR1, BCL2, and GATA3 were best predictors. RESULTS—All FFPE core biopsies gave interpretable data. Three-gene test predicted survival with high significance (p<2E-07).




## GENE EXPRESSION PROFILING FROM ARCHIVED FROZEN AND FIXED PROSTATE TISSUE AND URINE (Arul Chinniyayan, Univ. Michigan)

Seven studies to study expression of genes involved in prostate cancer progression and diagnosis. STUDY DESIGN— •Studies use WTA cDNA from frozen and fixed tissue, urine. •Array studies employ as few as 2,000 cells. •cDNA and oligonucleotide arrays used. •Expression results correlated to gene fusion and amplification events studied by WGA of the same samples. RESULTS—Genes and pathways most correlated with prostate cancer identified, as well as identification of urine expression biomarkers.




## THE CANCER GENOME ATLAS STUDY OF AMPLIFIED DNA ON CGH ARRAYS

Analysis of genomic copy number alterations using Agilent oligo 244K microarray and WGA. Below is a comparison of two aCGH profiles of the same GBM sample. Top profile has been obtained by a standard protocol required 2 mg gDNA. Lower profile is the result of GenomePlex WGA from 10 ng gDNA. Data from Dr. Alexei Protogopov of the Dana Farber Cancer Institute. (Published Sept. 2007 in TCGA Newsletter)



## SANGER INSTITUTE CGH STUDY OF SINGLE COPIES OF SORTED CHROMOSOMES


Single human chromosomes were sorted into wells of microplates. The wells were amplified at Rubicon using a single-cell GenomePlex kit. Amplified DNA from single wells were labeled and hybridized to BAC microarrays and as FISH probes. Translocations were mapped on the microarrays and chromosome spreads. Data from Dr. Nigel Carter.



Refs from Nigel Carter (Sanger Institute) 1. Gribble et al. (2004) Applications of Combined DNA Microarray and Chromosome Sorting Technologies. Chromosome Research 12:35-43. 2. Gribble et al. (2004) Chromosome Paints from Single Copy Chromosomes. Chromosome Research 12: 143-151. 3. Fieger et al. High Resolution Array CGH of Single Cells. (2006) NAR.

## CHROMATIN IMMUNOPRECIPITATION STUDY OF THE FIDELITY OF GENOMEPLEX AMPLIFICATION AND LINKER-MEDIATED PCR (Farnham lab)

Chromatin immunoprecipitation (ChIP) has proven to be a powerful tool, allowing the detection of protein-DNA interactions in living cells. A single reaction does not yield enough DNA to perform genomic profiling of protein binding, so amplification is necessary. GenomePlex was found to accurately represent the ChIP DNA much better than LM-PCR (Paggy Farnham laboratory)



Method	Total Peaks	Overlapping	Overlap (%)
LM-PCR	543	82	15
Post-PCR	491	343	70
WGA	443	280	63

- O'Geen et al. (2006) Comparison of Sample Preparation Methods for ChIP-chip Assays. BioTechniques 41:577-580.
- Acevedo, et al. (2007) Genome-Scale ChIP-chip Analysis Using 10,000 Human Cells. BioTechniques 43:791-797.

## MethylPlex and Bisulfite Test Concordance

MethylPlex Q-PCR assays are highly concordant with MethylLight assays for the same promoters. In this example of the promoter for APC1, the methylation indexes of the bisulfite-based and 2 MethylPlex assays were compared for 20 normal and twenty prostate cancer tissue samples. The analytical and clinical performance of both types of assays were very similar.



## METHYLPLEX CAN BE USED TO PROFILE METHYLATION ACROSS THE ENTIRE GENOME USING AGILENT OR NIMBLEGEN PROMOTER ARRAYS.

This slide shows ratio of methylation in a cancer tissue divided by the methylation of a normal tissue on an oligonucleotide array region with 5 b tiling. The differential methylation is detected across one of the CpG islands.



## METHYLPLEX CAN BE USED TO DISCOVER BIOMARKERS THAT CAN DISTINGUISH TWO SUBTYPES OF THIS CANCER FROM NORMAL, ADJACENT TISSUE FROM THE SAME PATIENTS.

This example is from a study of FFPE sections from the tumors and adjacent normal tissue. The plot shows two distinct clusters for cancer subtype 1 and cancer subtype 2, separate from normals.



## ILLUMINA AND GSK STUDY OF GENOTYPING ACCURACY OF GENOMEPLEX WGA

Genomic and WGA DNA were genotyped using the Illumina BeadArray™. Ten samples were analyzed at 2,320 sites. All data were pooled as there was no significant difference in number of calls or concordance among samples. The high concordance between gDNA and WGA DNA results demonstrate the accuracy and reproducibility of GenomePlex.

Assays Run	Genomic DNA	Amplified DNA
Calls Made	23,200	23,200
% Total	100%	100%
% Genomic Concordance*	99.8%	99.8%

- Barker et al. (2004) Two Methods of Whole-Genome Amplification Enable Accurate Genotyping Across a 2320 SNP Linkage Panel. Genome Research 901-907.
- Barnes et al., (2006) Polymorphisms in the novel gene acylxacyl hydrogenase (AOAH) are associated with asthma and associated phenotypes. J. Allergy and Clinical Immunology 118:70-77.
- Barnes et al. (2006) Variants in the gene encoding C3 are associated with asthma and related phenotypes among African Caribbean families. Genes and Immunity 7:27-35.

## USE OF GENOMEPLEX FOR ARRAY CGH OF MICRODISSECTED FFPE SECTIONS (Du lab, Cambridge UK)

Johnson et al. (2006) Application of Array CGH on Archival Formalin-Fixed Paraffin-Embedded Tissues Including Small Numbers of Microdissected Cells. Laboratory Investigation (2006) 1 - 11.

Sample	Maximum DNA Fragment (ng) amplified	Quality of aCGH by inspection	Pearson's correlation to aCGH data from corresponding frozen tissue	400 ng DNA from FFPE tissue	20 ng DNA from FFPE tissue with WGA	2000 cells from FFPE tissue with WGA
Xenografts						
X1	400	Good	0.98	0.96	0.91	
X2	400	Good	0.91	0.90	0.91	
X3	400	Good	0.88	0.88	0.91	
X4	400	Good	0.96	0.89	0.91	
X5	300	Good	0.83	0.82		
X6	100	Poor	0.72			
X7	100	Poor	0.77			
X8	100	Failed to label	0.50			
X9	100	Failed to label	0.51			
X10	100	Failed to label				
X11	100	Failed to label				
Primary glioblastoma						
P1	400	Good	0.91			
P2	400	Good	0.87			0.91
P3	300	Good	0.86			
P4	300	Good	0.91			