

ISBER Biospecimen Science Working Group

Chair : Fay Betsou

Co-chair : Elaine Gunter

2009-2010 members: G. Ashton, M. Barnes, E. Benson, D. Coppola, Y. DeSouza, J. Eliason, B. Glazer, F. Guadagni, D. Horsfall, C. Kleeberger, S. Lehmann, A. Prasad, K. Shea, A. Skubitz, S. Somiari


The Goal of the Working Group

- Biospecimen Science
 - Scientific knowledge about the biomolecular quality of samples
- Biospecimen Research

1st action

- **Biospecimen Research literature compilation**
 - Justification
 - A lot of data already published, but not always easy to find
 - Results :
 - **385 references (2009)**
 - **680 references (2010)**


Biospecimen Science Literature Compilation



<http://www.isber.org/wg/BS-WG-LitComp.html>


Human Biospecimen Science Studies in a clinical setting

Human Biospecimen Science Studies in a clinical setting




Human Biospecimen Science Studies in a research setting

Human Biospecimen Science Studies in a research setting



Environmental Biospecimen Science

Environmental Biospecimen Science



2nd action

- **Standard Biospecimen Research Experimental Protocol**

- Justification

- Extreme variability in experimental designs makes comparing results very difficult
- Example

- Results : experimental protocol including standardized options, for application in:

- Preclinical exploratory phase
- Clinical assay validation phase

Standard Biospecimen Research Experimental Protocol options

■ Plasma sample

Pre-centrifugation:

- Type of tube: **EDTA**, Heparin, ACD, Citrate, Sodium Fluoride/ Sodium Oxalate, Micro Tube with gel
- Delay time 1 (Δt_1): **<2h**, 2h, 4h, 8h, 24h, 48h, 72h
- Temperature: 20-25°C, **3-7°C**, 35-38°C

Centrifugation:

- Temperature: **20-22°C**, 3-7°C
- Number of centrifugations: 1, **2**

Post-centrifugation:

- Delay time 2 (Δt_2): **<1h**, 2h, 4h, 8h, 24h, 48h, 72h
- Pre-aliquoting temperature: 22-25°C, **3-7°C**, 35-38°C
- Storage temperature: -20°C, **-80°C**, vapor-phase liquid nitrogen
- Number of freeze-thaws: **1**, 2, 4, 10, 20, 30
- Storage duration (years): **1.5**, 3, 6, 9, 12, n
- Secondary container: **Polypropylene tube**, plastic straw, glass vial

Standard Biospecimen Research Experimental Protocol options

■ Urine sample

Pre-centrifugation:

- Type of collection: First void, Later void, 12h, **24h**
- Type of container: With protease inhibitors, **Without protease inhibitors**
- Delay time 1 (Δt_1): **<2h**, 2h, 4h, 8h, 24h, 48h
- Temperature: 20-25°C, **3-7°C**, 35-38°C, -80°C

Centrifugation:

- Temperature: **20-22°C**, 3-7°C
- g, Delay time 2 (Δt_2) : 1000g 10min, **2000g 10min**, 4000g 10min
- Number of centrifugations: 1, **2**

Post-centrifugation:

- Secondary container: **Polypropylene microtube**, Plastic straw, Polypropylene 5ml tube
- Delay time 3 (Δt_3): **<2h**, 2h, 4h, 8h, 24h, 48h, 72h
- Pre-aliquoting temperature: 22-25°C, **3-7°C**, 35-38°C
- Storage temperature: -20°C, **-80°C**, vapor-phase liquid nitrogen
- Number of freeze-thaws: **1**, 2, 4, 10, 20, 30
- Storage duration (years): **1**, 3, 6, 9, 12, n

Standard Biospecimen Research Experimental Protocol options

■ Solid tissue sample

Pre-fixation:

- Warm ischemia time: 0, **≤10min**, 15-20min, 25-30min, 50-60min
- Cold ischemia time (delay time 1 (Δt_1) before fixation): 2min, 15min, **30min**, 60min, 120min

Fixation:

- Type of fixation :
 - Cryopreservation: **Snap Freezing**, OCT embedding
 - No cryopreservation: RNA-Later, **Buffered formalin**, Zinc

Post-fixation:

- Delay time 2 (Δt_2) from OCT embedding to freezing: 0, 15min, **4-8h**, 24h, 48h, 72h
- Storage temperature (for cryopreservation): **-80°C**, vapor-phase liquid nitrogen
- Storage duration (years): **1**, 3, 6, 9, 12, n

Standard Biospecimen Research Experimental Protocol

- **Impact**

- Harmonization of Biospecimen research and comparability of results
- Identification and validation of new quality control tools
- Development of standard nomenclature for sample collection / processing / storage procedures

3rd action

- **Logistics and Sample Transport**

- Justification

- While *in-vitro* preanalytical factors can be controlled, shipping / logistics factors are out of the repository's control and may have severe impacts on the sample's quality

- Results

- Overview, recommendations, validation protocol

<http://www.isber.org/wg/BS-WG-LogTransp.html>

4th action

- **Standard PREanalytical Coding for Biospecimens : SPREC**
 - Justification
 - Traceability of biospecimen preanalytical variations
 - Results
 - Details on pre-analytical sample processing
 - 7-element-long biospecimen characterization code
 - Easy to implement
 - Flexible

SPREC-01

Fluid samples

TYPE OF SAMPLE	TYPE OF PRIMARY CONTAINER	PRE-CENTRIFUGATION	CENTRIFUGATION	SECOND CENTRIFUGATION	POST-CENTRIFUGATION	STORAGE
• ASC	• ACD	• RT <2h A	• RT 10mn <3000g nb A	• RT 10mn <3000g A	• <1h 3-7°C A	• PP 2ml -80°C A
• AMN	• CPD	• 3-7°C <2h B	• RT 10mn <3000g b B	• RT 10mn <3000g b B	• <1h RT B	• PP 2ml -20°C B
• BAL	• HEP	• RT 2-4h C	• 4°C 10mn <3000g C	• 4°C 10mn <3000g C	• 1-2h 3-7°C C	• Cryo 2ml LN C
• BLD	• HIR	• 3-7°C 2-4h D	• 4°C 10mn <3000g b D	• 4°C 10mn <3000g b D	• 1-2h RT D	• Cryo 2ml -80°C D
• BMA	• ORG	• RT 4-8h E	• RT 10mn 3-6000g b E	• RT 10mn 3-6000g b E	• 2-8h 3-7°C E	• Cryo 2ml PRF -135C E
• BMK	• PAX	• 3-7°C 4-8h F	• 4°C 10mn 3-6000g b F	• 4°C 10mn 3-6000g b F	• 2-8h RT F	• Straw LN F
• BUC	• PED	• RT 8-12h G	• RT 10mn 6-10000g b G	• RT 10mn 6-10000g b G	• 8-24h 3-7°C G	• Straw -80°C G
• BUF	• PI1	• 3-7°C 8-12h H	• 4C 10mn 6-10000g b H	• 4°C 10mn 6-10000gbH	• 8-24h RT H	• Straw -20°C H
• CEL	• PIX	• RT 12-24h I	• RT 10mn >10000g b I	• RT 10mn >10000g b I	• >24h 3-7°C I	• Straw PRF -135°C I
• CEN	• PPS	• 3-7°C 12-24h J	• 4°C 10mn >10000g b J	• 4°C 10mn >10000g b J	• >24h RT J	• PP 5ml -80°C J
• CLN	• PXD	• RT 24-48h K	• No centrifugation N	• No 2 nd centrifugation	• X	• PP 5ml -20°C K
• CRD	• PXR	• 3-7°C 24-48h L	• X	N	• Z	• Mplate -80°C L
• CSF	• SCI	• RT >48h M	• Z	• X		• Mplate -20°C M
• NAS	• SED	• 3-7°C >48h N		• Z		• Paraffin RT P
• PEL	• SPO	• 35-38°C <2h O				• X
• PEN	• SST					• Z
• PFL	• TEM					
• PL1	• TRC					
• PL2	• XXX					
• SAL	• ZZZ					
• SEM						
• SER						
• SPT						
• STL						
• SYN						
• TER						
• U24						
• URN						
• ZZZ						

SPREC-01

Solid samples

TYPE OF SAMPLE	TYPE OF COLLECTION	WARM ISCHEMIA TIME	COLD ISCHEMIA TIME	FIXATION TYPE	FIXATION TIME	STORAGE
<ul style="list-style-type: none"> •CEN •CLN •FNA •HAR •LCM •PEN •TIS •LCM •ZZZ 	<ul style="list-style-type: none"> •A06 •A12 •A24 •A48 •A72 •BPS •FNA •PUN •SRG •SWB •ZZZ 	<ul style="list-style-type: none"> •<2min A •2-10min B •10-20min C •20-30min D •30-60min E •>60min F •X •N •Z 	<ul style="list-style-type: none"> •<2min A •2-10min B •10-20min C •20-30min D •30-60min E •>60min F •X •N •Z 	<ul style="list-style-type: none"> •ACA •ALD •ETH •FOR •SNP •NAA •NBF •OCT •RNL •XXX •ZZZ 	<ul style="list-style-type: none"> •<15min A •15min-1h B •1-4h C •4-8h D •8-24h E •24-48h F •48-72h G •X •Z 	<ul style="list-style-type: none"> •PP 2ml -80°C A •PP 2ml -20°C B •Cryo 2ml LN C •Cryo 2ml -80°C D •Cryo 2ml PRF -135°C E •Straw LN F •Straw -80°C G •Straw -20°C H •Straw PRF -135°C I •PP 5ml -80°C J •PP 5ml -20°C K •Mplate -80°C L •Mplate -20°C M •Paraffin RT P •X •Z


SPREC-01

Scope of application

- Primary samples
 - Specimens directly collected from the donor
- Simple derivatives
 - Samples prepared by simple laboratory manipulation
- Complex derivatives, out of scope
 - Cell disruption
 - Cell selection
 - Multi-step chemical manipulation
 - Acidification
 - Digestion
 - Precipitation
 - Deproteinisation
 - Desalting

SPREC-01

Biospecimen description examples



SPREC-01

Impact

- Get incorporated in databases and QMS
- Facilitate and consolidate biomarker identification research
- Support biospecimen research
- Support laboratory accreditation
- Complement reporting recommendations
 - STARD
 - « baseline data »
 - « baseline demographic and clinical characteristics »
 - www.consort-statement.org/initiatives/newstand.htm
 - STROBE
 - « descriptive data »
 - www.strobe-statement.org/index.html
 - REMARK
 - « type of biological material used, methods of preservation and storage »
 - [J Nat Cancer Inst 2005;97:1180-84](#)

Future actions

- Proficiency Testing programs development for assessment of Quality Control assay performance
- Inter-laboratory evaluation of novel technologies for sample collection, processing, storage, transport

Thank you!!!

ISBER Biospecimen Science Working Group



Amy Skubitz



David Horsfall



Elaine Gunter



Kathi Shea



Mike Barnes



Cynthia Kleeberger




James Eliason



Stella Somiari



Erica Benson



Sylvain Lehmann



Yvonne DeSouza



Domenico Coppola



Barbara Glazer



Fay Betsou



Fiorella Guadagni



Anil Prasad



Garry Ashton