



Glycoproteomic analysis of formalin-fixed and paraffin-embedded tissues and OCT-embedded frozen tissues

Yuan Tian, Kay Gurley, Steve G. Bova, Christopher J. Kemp, and Hui Zhang

Introduction

- ❖ FFPE tissues and OCT-embedded tissues are commonly used and archived with pathological and clinical information;
- ❖ Applications:
 - Histology
 - Immunohistochemistry for molecular marker
 - In situ hybridization
 - Proteomics

FFPE: Formalin- Fixed Paraffin-Embedded

OCT: Optimal Cutting Temperature compound

Introduction

Glycoproteins as Serum Tumor Markers

- Prevalent in extracellular proteins
- Most likely detected as surrogate markers
- Altered glycosylation associated with cancer
- All FDA approved protein serum markers are glycoproteins

Introduction

❖ Can these processed tissues be used for proteomics analysis by mass spectrometer?

OCT compound consists of a resin-polyvinyl alcohol, benzalonium chloride, an antifungal agent, and polyethylene glycol to lower the freezing temperature


❖ Glycosylation is preserved in the embedded tissues?

Aim

- To investigate the feasibility of glycoproteomic analysis using FFPE tissue and OCT-embedded frozen tissue

Method

Solid phase extraction of glycopeptides (SPEG) and quantitative glycoproteomics




Zhang, H. Li, X-J, Martin, D. M., and Aebersold, R. *Nature Biotechnology* (2003).

Tian, Y., Zhou, Y., Elliott, S., Aebersold, R., and Zhang, H. *Nature Protocols* (2007).


Zhou, Y. et al. Isolation of *N*-linked glycopeptides from plasma. *Analytical Chemistry* (2007).

Zou, Z., Ibisate, M., Zhou, Y., Aebersold, R. Xia, Y., and Zhang, H. Synthesis and Evaluation of Superparamagnetic Nanoparticles for Extraction of Glycopeptides in the Microtiter Plate Format *Analytical Chemistry* (2008).



Breast cancer tissues/cells
1208 glycosites

B and T lymphocytes
407 glycosites




Zhang, H. et al. UniPep, a database for human N-linked glycosites: A Resource for Biomarker Discovery. *Genome Biology* (2006).

Zhang, H. et al. Mass spectrometric detection of tissue-derived proteins in blood. *Molecular & Cellular Proteomics* (2007)

Studies


I. Glycoprotein analysis of FFPE tissue vs. frozen tissue



Studies

I. Glycoprotein analysis of FFPE tissue vs. frozen tissue

1. Unique glycosylation sites



Tian Y. et al *J Proteome Res* 2009, 8 (4), 1657-1662

Studies

I. Glycoprotein analysis of FFPE tissue vs. frozen tissue

2. Tryptic ends & missed cleavage sites


	% of 1 tryptic end	% of 2 tryptic end	% of -R	% of -K	missed cleavage site %
FFPE lung	63.10	36.9	51.63	48.37	16.07
Frozen lung	64.68	35.32	50.43	49.57	17.06

Tian Y. et al *J Proteome Res* 2009, 8 (4), 1657-1662

Studies

I. Glycoprotein analysis of FFPE tissue vs. frozen tissue

3. Cellular location of identified glycoproteins




Tian Y. et al *J Proteome Res* 2009, 8 (4), 1657-1662

Studies

II. Identification of glycoproteins associated with aggressive prostate cancer using OCT-embedded tissues


- Prostate cancer is the most common malignancy in men in the United States;
- No good method to reliably distinguish aggressive from non-aggressive prostate cancer




Studies

II. Identification of glycoproteins associated with aggressive prostate cancer using OCT-embedded tissues

A
 microfibrillar-associated protein 4_pep1




B
 periostin




Studies

II. Identification of glycoproteins associated with aggressive prostate cancer using OCT-embedded tissues

C




D



Studies

II. Identification of glycoproteins associated with aggressive prostate cancer using OCT-embedded tissues

Verification using Western blot



Studies

II. Identification of glycoproteins associated with aggressive prostate cancer using OCT-embedded tissues

Verification using immunohistochemistry




Cathepsin L

H&E staining




100X

600X

Normal



Primary prostate cancer




Studies

II. Identification of glycoproteins associated with aggressive prostate cancer using OCT-embedded tissues


Verification using immunohistochemistry

Cathepsin L


LN Met



Liver Met



Pericardial Bulk
Met




Spleen Met



Posterior Subdura
Met



Normal



Studies

II. Identification of glycoproteins associated with aggressive prostate cancer using OCT-embedded tissues

Verification using immunohistochemistry




Periostin

H&E




100X

600X

Normal



Primary prostate cancer




Studies

II. Identification of glycoproteins associated with aggressive prostate cancer using OCT-embedded tissues


Verification using immunohistochemistry

Periostin


LN Met




Liver Met



Prostate cancer



Normal Prostate



Conclusion

- ❖ FFPE tissues and OCT-embedded frozen tissues are feasible to be used in glycoproteomics analysis
- ❖ Glycans can be used as tags to isolate glycoproteins from FFPE and OCT-embedded biospecimens

Acknowledgements

Johns Hopkins University

Hui Zhang

Daniel Chan

Xiaer Sun

G. Steve Bova

Dylan Bach

Robert N. Cole

Robert O'Meally

Fred Hutchinson Cancer Research Center

Christopher J. Kemp

Kay Gurley

Grant support

NIH/NCI

Early Detection Research

Network (EDRN)

Patrick C. Walsh Prostate Cancer
Research Fund