

**Biospecimen Pre-Analytical Variables
 (BPV) Formalin-Fixed, Paraffin-Embedded
 (FFPE) Section/Stain Form**

OP-0004-F3

VER. 03.01

Effective Date: 11/10/2014

Page 1 of 5

BPV Case ID _____

Affix BPV Case ID Barcode Label

Form Completed By: _____

Date Form Was Completed: __/__/____
 (MM/DD/YYYY)

Biospecimen Source Site:

Primary Organ: _____

Tissue Bank ID: _____

FFPE Sectioning

Slide Prep Technician Name:

Slides Created

Module

Slide ID:

Slide ID:

Slide ID:

Slide ID:

Slide ID:

Slide ID:

Slide ID:

Slide ID:

Biospecimen Pre-Analytical Variables (BPV) Formalin-Fixed, Paraffin-Embedded (FFPE) Section/Stain Form

OP-0004-F3

VER. 03.01

Effective Date: 11/10/2014

Page 2 of 5

BPV Case ID: _____

Tissue Bank ID: _____

Form Completed By: _____

Date Completed ____/____/____
MM/DD/YYYY

FFPE Sectioning Details

1. Slide Prep Standard of Protocol (SOP):		Your SOP governing preparation of FFPE tissue sections
2. Microtome:	Select one: <input type="radio"/> Microm Ergostar HM200 <input type="radio"/> Other, specify:	
3. Microtome Blade Type:	Select one: <input type="radio"/> Low-profile Sakura Accu-Edge <input type="radio"/> Other, specify:	
4. Microtome Blade Age:	Select one: <input type="radio"/> Fresh <input type="radio"/> Other, specify:	
5. Preparation of Block Face for Sectioning:	Select one: <input type="radio"/> 20 minutes on ice <input type="radio"/> Other, specify:	
6. Section Thickness:	Select one: <input type="radio"/> 4–5 micrometers <input type="radio"/> Other; specify:	

**Biospecimen Pre-Analytical Variables
 (BPV) Formalin-Fixed, Paraffin-Embedded
 (FFPE) Section/Stain Form**

OP-0004-F3

VER. 03.01

Effective Date: 11/10/2014

Page 3 of 5

BPV Case ID: _____

Tissue Bank ID: _____

Form Completed By: _____

Date Completed ____/____/____
MM/DD/YYYY

7. Slide Charge:	Select one: <input type="radio"/> Charged <input type="radio"/> Uncharged <input type="radio"/> Other, specify:	
8. Water Bath Temperature:	Select one: <input type="radio"/> 41–45°C <input type="radio"/> Other, specify:	
9. Microtome Maintenance:	Select one: <input type="radio"/> Daily <input type="radio"/> Other, specify:	
10. Waterbath Maintenance:	Select one: <input type="radio"/> Daily <input type="radio"/> Other, specify:	
11. Any additional comments related to preparation of FFPE tissue sections:		

**Biospecimen Pre-Analytical Variables
 (BPV) Formalin-Fixed, Paraffin-Embedded
 (FFPE) Section/Stain Form**

OP-0004-F3

VER. 03.01

Effective Date: 11/10/2014

Page 4 of 5

BPV Case ID: _____

Tissue Bank ID: _____

Form Completed By: _____

Date Completed ____/____/____
 MM/DD/YYYY

Hematoxylin and Eosin (H&E) Staining

<p>12. H&E Stain SOP:</p>		<p>Your SOP governing H&E staining and slide mounting of tissue sections</p>
<p>13. H&E Time in Oven:</p>	<p>Select one: <input type="radio"/> 20 minutes <input type="radio"/> Other, specify:</p>	
<p>14. H&E Oven Temperature:</p>	<p>Select one: <input type="radio"/> 60°C <input type="radio"/> Other, specify:</p>	
<p>15. H&E De-Paraffin Method:</p>	<p>Select one: <input type="radio"/> Manual <input type="radio"/> Automated stainer <input type="radio"/> Other, specify:</p>	
<p>16. H&E Stain Method:</p>	<p>Select one: <input type="radio"/> Manual <input type="radio"/> Automated stainer <input type="radio"/> Other, specify:</p>	

**Biospecimen Pre-Analytical Variables
 (BPV) Formalin-Fixed, Paraffin-Embedded
 (FFPE) Section/Stain Form**

OP-0004-F3

VER. 03.01

Effective Date: 11/10/2014

Page 5 of 5

BPV Case ID: _____

Tissue Bank ID: _____

Form Completed By: _____

Date Completed ____/____/_____
MM/DD/YYYY

17. H&E Clearing Method:

Select one:
 Manual
 Automated stainer
 Other, specify:

18. H&E Cover Slipping:

Select one:
 Manual
 Other, specify:

19. H&E Equipment Maintenance:

Select one:
 Daily, weekly, bimonthly, or per SOP
 Other, specify:

20. Additional comments related to preparation of FFPE H&E-stained slides: